

USO1 Rabbit mAb

Catalog No: #57097

Package Size: #57097-1 50ul #57097-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	USO1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF
Species Reactivity	Human Mouse Rat
Specificity	USO1 Antibody detects endogenous levels of total USO1
Immunogen Description	A synthesized peptide derived from human USO1
Other Names	P115; TAP; Uso1; VDP;
Accession No.	Uniprot:O60763
Uniprot	O60763
Calculated MW	108kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB: 1:500~1:2000 IHC: 1:50~1:200 ICC/IF: 1:50~1:200

Images



Western blot analysis of USO1 expression in HepG2 cell lysate.

Product Description

General vesicular transport factor required for intercisternal transport in the Golgi stack; it is required for transcytotic fusion and/or subsequent binding of the vesicles to the target membrane.

Note: This product is for in vitro research use only