

## NUDT19 Rabbit mAb

Catalog No: #57098

Package Size: #57098-1 50ul #57098-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	NUDT19 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF
Species Reactivity	Human
Specificity	NUDT19 Antibody detects endogenous levels of total NUDT19
Immunogen Description	A synthesized peptide derived from human NUDT19
Other Names	NUDT19; RP2;
Accession No.	Uniprot:A8MXV4
Uniprot	A8MXV4
Calculated MW	42kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

## Application Details

WB: 1:500~1:2000 IHC: 1:50~1:200 ICC/IF: 1:50~1:200

## Images



Western blot analysis of NUDT19 expression in 293T cell lysate.

## Product Description

Coenzyme A diphosphatase that mediates the hydrolysis of a wide range of CoA esters, including choloyl-CoA and branched-chain fatty-acyl-CoA esters. At low substrate concentrations medium and long-chain fatty-acyl-CoA esters are the primary substrates.

---

Note: This product is for in vitro research use only