

TOR4A Rabbit mAb

Catalog No: #57099

Package Size: #57099-1 50ul #57099-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	TOR4A Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human
Specificity	TOR4A Antibody detects endogenous levels of total TOR4A
Immunogen Description	A synthesized peptide derived from human TOR4A
Other Names	C9orf167;
Accession No.	Uniprot:Q9NXH8
Uniprot	Q9NXH8
Calculated MW	47kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB: 1:500~1:2000

Images



Western blot analysis of TOR4A expression in HeLa cell lysate.

Product Description

ATP-binding.

Note: This product is for in vitro research use only