

MST4 Rabbit mAb

Catalog No: #57103

Package Size: #57103-1 50ul #57103-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	MST4 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human Mouse Rat
Specificity	MST4 Antibody detects endogenous levels of total MST4
Immunogen Description	A synthesized peptide derived from human MST4
Other Names	MASK; Mst4; STE20 like kinase MST4; STK26;
Accession No.	Uniprot:Q9P289
Uniprot	Q9P289
Calculated MW	50kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB: 1:500~1:2000 IHC: 1:50~1:200

Images



Western blot analysis of MST4 expression in HeLa cell lysate.

Product Description

Mediator of cell growth. Modulates apoptosis.

Note: This product is for in vitro research use only