

GALT Rabbit mAb

Catalog No: #57104

Package Size: #57104-1 50ul #57104-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description


Product Name	GALT Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF
Species Reactivity	Human Mouse Rat
Specificity	GALT Antibody detects endogenous levels of total GALT
Immunogen Description	A synthesized peptide derived from human GALT
Other Names	GALT;
Accession No.	Uniprot:P07902
Uniprot	P07902
Calculated MW	43kDa
Concentration	0.2mg/ml
Formulation	Rabbit IgG in 10mM phosphate buffered saline , pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000

ICC/IF:1:50~1:200

Images



Western blot analysis of GALT expression in K562 cell lysate.

Product Description

Plays an important role in galactose metabolism.

Note: This product is for in vitro research use only