

DAZAP1 Rabbit mAb

Catalog No: #57107

Package Size: #57107-1 50ul #57107-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	DAZAP1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	DAZAP1 Antibody detects endogenous levels of total DAZAP1
Immunogen Description	A synthesized peptide derived from human DAZAP1
Other Names	DAZ associated protein 1; Dazap1;
Accession No.	Uniprot:Q96EP5
Uniprot	Q96EP5
Calculated MW	43kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB: 1:500~1:2000

Images



Western blot analysis of DAZAP1 expression in Jurkat cell lysate.

Product Description

RNA-binding protein, which may be required during spermatogenesis.

Note: This product is for in vitro research use only