

PPP2R5D Rabbit mAb

Catalog No: #57112

Package Size: #57112-1 50ul #57112-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	PPP2R5D Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP FC
Species Reactivity	Human Mouse Rat
Specificity	PPP2R5D Antibody detects endogenous levels of total PPP2R5D
Immunogen Description	A synthesized peptide derived from human PPP2R5D
Other Names	B56D; PPP2R5D; Tex271;
Accession No.	Uniprot:Q14738
Uniprot	Q14738
Calculated MW	60,65kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IHC:1:50~1:200ICC/IF:1:50~1:200IP:1:50FC:1:50

Images



Western blot analysis of PPP2R5D expression in HeLa cell lysate.

Product Description

The B regulatory subunit might modulate substrate selectivity and catalytic activity, and also might direct the localization of the catalytic enzyme to a particular subcellular compartment.

Note: This product is for in vitro research use only