

SNX5 Rabbit mAb

Catalog No: #57114

Package Size: #57114-1 50ul #57114-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	SNX5 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF
Species Reactivity	Human Mouse Rat
Specificity	SNX5 Antibody detects endogenous levels of total SNX5
Immunogen Description	A synthesized peptide derived from human SNX5
Other Names	SNX5; Sorting nexin 5;
Accession No.	Uniprot:Q9Y5X3
Uniprot	Q9Y5X3
Calculated MW	47kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IHC:1:50~1:200ICC/IF:1:50~1:200

Images



Western blot analysis of SNX5 expression in Jurkat cell lysate.

Product Description

May be involved in several stages of intracellular trafficking. Plays a role in macropinocytosis. Plays a role in the internalization of EGFR after EGF stimulation.

Note: This product is for in vitro research use only