

PGM1 Rabbit mAb

Catalog No: #57116

Package Size: #57116-1 50ul #57116-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	PGM1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF IP
Species Reactivity	Human Mouse Rat
Specificity	PGM1 Antibody detects endogenous levels of total PGM1
Immunogen Description	A synthesized peptide derived from human PGM1
Other Names	CDG1T; Glucose phosphomutase 1; GSD14; PGM1; Phosphoglucomutase 1;
Accession No.	Uniprot:P36871
Uniprot	P36871
Calculated MW	61kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000ICC/IF:1:50~1:200IP:1:50

Images



Western blot analysis of PGM1 expression in (1) 293 cell lysate; (2) RAW 264.7 cell lysate.

Product Description

This enzyme participates in both the breakdown and synthesis of glucose.

Note: This product is for in vitro research use only