

TNFAIP8 Rabbit mAb

Catalog No: #57117

Package Size: #57117-1 50ul #57117-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TNFAIP8 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human
Specificity	TNFAIP8 Antibody detects endogenous levels of total TNFAIP8
Immunogen Description	A synthesized peptide derived from human TNFAIP8
Other Names	NDED; SCCS2; TNFAIP8;
Accession No.	Uniprot:O95379
Uniprot	O95379
Calculated MW	18,21kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4ℳ C short term. Store at -20ℳ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000

Images



Western blot analysis of TNFAIP8 expression in A549 cell lysate.

Product Description

Acts as a negative mediator of apoptosis and may play a role in tumor progression. Suppresses the TNF-mediated apoptosis by inhibiting caspase-8 activity but not the processing of procaspase-8, subsequently resulting in inhibition of BID cleavage and caspase-3 activation.

Note: This product is for in vitro research use only