

## NOLC1 Rabbit mAb

Catalog No: #57120

Package Size: #57120-1 50ul #57120-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	NOLC1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF
Species Reactivity	Human
Specificity	NOLC1 Antibody detects endogenous levels of total NOLC1
Immunogen Description	A synthesized peptide derived from human NOLC1
Other Names	NOLC1; NOPP130; Nopp140; NS5ATP13; Nucleolar phosphoprotein p130; P130;
Accession No.	Uniprot:Q14978
Uniprot	Q14978
Calculated MW	110kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 $\Lambda$ C short term. Store at -20 $\Lambda$ C long term. Avoid freeze / thaw cycle.

## Application Details

WB:1:500~1:2000IHC:1:50~1:200ICC/IF:1:50~1:200

## Images



Western blot analysis of NOLC1 expression in HeLa cell lysate.

## Product Description

Related to nucleogenesis, may play a role in the maintenance of the fundamental structure of the fibrillar center and dense fibrillar component in the nucleolus. It has intrinsic GTPase and ATPase activities. May play an important role in transcription catalyzed by RNA polymerase I.

---

Note: This product is for in vitro research use only