

EMC4 Rabbit mAb

Catalog No: #57121

Package Size: #57121-1 50ul #57121-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	EMC4 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC IP
Species Reactivity	Human Mouse Rat
Specificity	EMC4 Antibody detects endogenous levels of total EMC4
Immunogen Description	A synthesized peptide derived from human EMC4
Other Names	EMC4; PIG17; tmem85; Transmembrane protein 85;
Accession No.	Uniprot:Q5J8M3
Uniprot	Q5J8M3
Calculated MW	18kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB: 1:500~1:2000 IHC: 1:50~1:200 IP: 1:50

Images



Western blot analysis of EMC4 expression in 293 cell lysate.

Product Description

May mediate anti-apoptotic activity.

Note: This product is for in vitro research use only