

DNAJA4 Rabbit mAb

Catalog No: #57123

Package Size: #57123-1 50ul #57123-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	DNAJA4 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IP
Species Reactivity	Human Mouse Rat
Specificity	DNAJA4 Antibody detects endogenous levels of total DNAJA4
Immunogen Description	A synthesized peptide derived from human DNAJA4
Other Names	DNAJA4; MST104; MSTP104; PRO1472;
Accession No.	Uniprot:Q8WW22
Uniprot	Q8WW22
Calculated MW	45kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IP:1:50

Images



Western blot analysis of DNAJA4 expression in HeLa cell lysate.

Product Description

Contains 1 CR-type zinc finger. Contains 1 J domain.

Note: This product is for in vitro research use only