

FAM50A Rabbit mAb

Catalog No: #57124

Package Size: #57124-1 50ul #57124-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	FAM50A Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF
Species Reactivity	Human Mouse Rat
Specificity	FAM50A Antibody detects endogenous levels of total FAM50A
Immunogen Description	A synthesized peptide derived from human FAM50A
Other Names	fam50a; HXC26; XAP5;
Accession No.	Uniprot:Q14320
Uniprot	Q14320
Calculated MW	40kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB: 1:500~1:2000 IHC: 1:50~1:200 ICC/IF: 1:50~1:200

Images



Western blot analysis of FAM50A expression in (1) HeLa cell lysate; (2) Raw264.7 cell lysate.

Product Description

May be a DNA-binding protein or transcriptional factor.

Note: This product is for in vitro research use only