

TAGAP Rabbit mAb

Catalog No: #57125

Package Size: #57125-1 50ul #57125-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	TAGAP Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IP
Species Reactivity	Human Mouse Rat
Specificity	TAGAP Antibody detects endogenous levels of total TAGAP
Immunogen Description	A synthesized peptide derived from human TAGAP
Other Names	ARHGAP47; IDDM21; Tagap; TAGAP1;
Accession No.	Uniprot:Q8N103
Uniprot	Q8N103
Calculated MW	65,90kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB: 1:500~1:2000 IP: 1:50

Images



Western blot analysis of TAGAP expression in Jurkat cell lysate.

Product Description

May function as a GTPase-activating protein and may play important roles during T-cell activation.

Note: This product is for in vitro research use only