

## FITC mouse IgG1,k Isotype Control

Catalog No: #IC002FITC

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	FITC mouse IgG1,k Isotype Control
Brief Description	FITC mouse IgG1,k Isotype Control
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1,k
Purification	Purified by Protein A.
Applications	FC IF
Concentration	10 <sup>3</sup> EOI/Test
Storage	Phosphate-buffered solution, pH 7.4, containing 0.09% sodium azide and 0.2% (w/v) BSA,Store at 4°C.DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.

## Background

Isotype controls are negative controls that are designed to measure the amount of non-specific background signal caused by the primary antibody. Isotype Control have no specificity for the target protein, yet retain all of the non-specific characteristics of an antibody.

Note: This product is for in vitro research use only