

TGR5 Polyclonal Antibody

Catalog No: #44221

Package Size: #44221-1 50ul #44221-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TGR5 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	Human; Mouse; Rat
Species Reactivity	WB;ELISA;IHC-p
Specificity	TGR5 Polyclonal Antibody detects endogenous levels of TGR5 protein.
Immunogen Description	The antiserum was produced against synthesized peptide derived from human GPBAR.
Accession No.	Swiss Prot:Q8TDU6GeneID:151306
Uniprot	Q8TDU6
Calculated MW	35kD
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

WB 1:500-2000;IHC-p 1:50-300; ELISA 2000-20000

Images



Immunohistochemical analysis of paraffin-embedded human tonsil.

Note: This product is for in vitro research use only