

Rabbit Anti-Human IgM

Catalog No: #L35658



Package Size: #L35658-1 100ul #L35658-2 500ul #L35658-3 1ml

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|-----------------------|---|
| Product Name | Rabbit Anti-Human IgM |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Purification | The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads. |
| Applications | E, WB, IH |
| Species Reactivity | Human |
| Immunogen Description | Human IgM |
| Formulation | Liquid in 0.01M Phosphate Buffered Saline, pH7.2, 0.05% Sodium Azide |
| Storage | Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles. |

Application Details

E (1/5000 - 1/20000), WB (1/5000 - 1/20000), IH (1/500 - 1/1000)

Note: This product is for in vitro research use only