

Human Ras-related protein Rab-26 (RAB26) ELISA Kit



Catalog No: #EK20477

Orders: order@signalwayantibody.com

Package Size: #EK20477-1 48T #EK20477-2 96T

Support: tech@signalwayantibody.com

Description

Product Name	Human Ras-related protein Rab-26 (RAB26) ELISA Kit
Applications	ELISA
Species Reactivity	Human

Application Details

Detect Range:125-8000 pg/mL

Sensitivity:51pg/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only