

## Thrombospondin 4 Polyclonal Antibody

Catalog No: #44224

Package Size: #44224-1 50ul #44224-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Thrombospondin 4 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB;IHC-p;IF/ICC;ELISA
Species Reactivity	Human; Mouse;Rat
Specificity	Thrombospondin 4 Polyclonal Antibody detects endogenous levels of Thrombospondin 4 protein.
Immunogen Description	The antiserum was produced against synthesized peptide derived from human THBS4.
Accession No.	Swiss Prot:P35443GeneID:7060
Calculated MW	106kD
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


## Application Details

WB 1:500-2000;


IHC-p 1:100-500;IF/ICC 1:100-500;

ELISA 1:5000-20000

## Images



Western blot analysis of 293T lysis using Thrombospondin 4 antibody.



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma, using THBS4 Antibody.

---

Note: This product is for in vitro research use only