

# Human Neuropeptide Y (NPY) ELISA Kit

Catalog No: #EK20486



Package Size: #EK20486-1 48T #EK20486-2 96T

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human Neuropeptide Y (NPY) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human

## Application Details

Detect Range:3.12-200pg/mL

Sensitivity:0.97pg/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 &mu;L

Assay Time:1-4.5h

Detection wavelength:450 nm

## Product Description

Detection Method:competitive inhibition

Note: This product is for in vitro research use only