

PHF5A rabbit pAb

Catalog No: #44277



Package Size: #44277-1 50ul #44277-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	PHF5A rabbit pAb
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	WB;ELISA;IHC
Species Reactivity	Human; Mouse;Rat
Specificity	This antibody detects endogenous levels of PHF5A at Human/Mouse/Rat
Immunogen Description	Synthesized peptide derived from human PHF5A
Accession No.	Swiss-Prot#: Q7RTV0 NCBI Gene#: 84844
Calculated MW	12kD
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

WB 1:500-2000

IHC-p 1:50-300

ELISA 2000-20000

Note: This product is for in vitro research use only