

Tbp7 Rabbit mAb

Catalog No: #57143

Package Size: #57143-1 50ul #57143-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Tbp7 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF
Species Reactivity	Human Mouse Rat
Specificity	Tbp7 Antibody detects endogenous levels of total Tbp7
Immunogen Description	A synthesized peptide derived from human Tbp7
Other Names	MIP224; PSMC4; RPT3; Tat binding protein 7; TBP 7;
Accession No.	P43686
Calculated MW	47kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

WB 1:500~1:2000 ICC/IF 1:50~1:200

Images



Western blot analysis of Tbp7 expression in HeLa cell lysate.

Background

The 26S protease is involved in the ATP-dependent degradation of ubiquitinated proteins. The regulatory (or ATPase) complex confers ATP dependency and substrate specificity to the 26S complex.

Note: This product is for in vitro research use only