

## hCG beta Rabbit mAb

Catalog No: #57163



Package Size: #57163-1 50ul #57163-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	hCG beta Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Purification	Affinity-chromatography
Applications	IHC
Species Reactivity	Human
Specificity	hCG beta Antibody detects endogenous levels of total hCG beta
Other Names	CG beta; CgB; CGB3; CGB5; CGB7; CGB8; Choriongonadotropin subunit beta 3; Choriongonadotropin subunit beta 7; Choriongonadotropin subunit beta; Chorionic gonadotropin beta polypeptide; Chorionic gonadotropin chain beta;
Accession No.	P0DN87/Q6NT52
Concentration	1 mg/ml
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

## Application Details

IHC 1:50~1:200

## Background

Function: Stimulates the ovaries to synthesize the steroids that are essential for the maintenance of pregnancy. Tissue specificity: Placenta. Similarity: Belongs to the glycoprotein hormones subunit beta family. Developmental stage: Made by the first trimester placenta.

Note: This product is for in vitro research use only