

IL13 Rabbit mAb

Catalog No: #57165

Package Size: #57165-1 50ul #57165-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	IL13 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human
Specificity	IL13 Antibody detects endogenous levels of total IL13
Immunogen Description	A synthesized peptide derived from human IL13
Other Names	ALRH; BHR1; Il13; Interleukin 13; NC30; P600;
Accession No.	P35225
Calculated MW	13kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

WB 1:500~1:2000

Images



Western blot analysis of IL13 expression in human recombinant protein cell lysate.

Background

Cytokine (PubMed:8096327, PubMed:8097324). Inhibits inflammatory cytokine production (PubMed:8096327). Synergizes with IL2 in regulating interferon-gamma synthesis (PubMed:8096327).

Note: This product is for in vitro research use only