

Paxillin (Phospho-Tyr118) Conjugated Antibody

Catalog No: #C14217



Package Size: #C14217-AF350 100ul #C14217-AF405 100ul #C14217-AF488 100ul
 #C14217-AF555 100ul #C14217-AF594 100ul #C14217-AF647 100ul
 #C14217-AF680 100ul #C14217-AF750 100ul #C14217-Biotin 100ul

Orders: order@signalwayantibody.com
 Support: tech@signalwayantibody.com

Description

Product Name	Paxillin (Phospho-Tyr118) Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Human Mouse Rat
Specificity	Phospho-Paxillin (Y118) Antibody detects endogenous levels of Phospho-Paxillin (Y118)
Immunogen Description	A synthesized peptide derived from human Paxillin
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	Paired box protein Pax 1; PAX 1; PAX1; Paxillin alpha; Paxillin; PXN; PXN protein;
Accession No.	Uniprot:P49023
Uniprot	P49023
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	68kDa
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Product Description

Cytoskeletal protein involved in actin-membrane attachment at sites of cell adhesion to the extracellular matrix (focal adhesion).

Note: This product is for in vitro research use only