

## Histone H3.3 (Phospho-Ser31) Conjugated Antibody

Catalog No: #C14268



Package Size: #C14268-AF350 100ul #C14268-AF405 100ul #C14268-AF488 100ul  
 #C14268-AF555 100ul #C14268-AF594 100ul #C14268-AF647 100ul  
 #C14268-AF680 100ul #C14268-AF750 100ul #C14268-Biotin 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
 Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Histone H3.3 (Phospho-Ser31) Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	WB IHC ICC/IF FC
Specificity	Phospho-Histone H3.3 (S31) Antibody detects endogenous levels of total Phospho-Histone H3.3 (S31)
Immunogen Description	A synthesized peptide derived from human Phospho-Histone H3.3 (S31)
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	Histone H3.3;
Accession No.	Uniprot:P84243
Uniprot	P84243
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	15kDa
Storage	Store at 4°C in dark for 6 months

## Application Details

## Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

## Product Description

---

Variant histone H3 which replaces conventional H3 in a wide range of nucleosomes in active genes. Constitutes the predominant form of histone H3 in non-dividing cells and is incorporated into chromatin independently of DNA synthesis.

---

Note: This product is for in vitro research use only