

# Nucleophosmin (Phospho-Thr199) Conjugated Antibody

Catalog No: #C14269

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #C14269-AF350 100ul #C14269-AF405 100ul #C14269-AF488 100ul

#C14269-AF555 100ul #C14269-AF594 100ul #C14269-AF647 100ul

#C14269-AF680 100ul #C14269-AF750 100ul #C14269-Biotin 100ul

## Description

Product Name	Nucleophosmin (Phospho-Thr199) Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	WB IHC ICC/IF FC
Specificity	Phospho-Nucleophosmin (T199) Antibody detects endogenous levels of total Phospho-Nucleophosmin (T199)
Immunogen Description	A synthesized peptide derived from human Phospho-Nucleophosmin (T199)
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	NO38; NPM; NPM1; Nucleolar phosphoprotein B23; Nucleophosmin; nucleophosmin nucleoplasmin family member 1; Numatrin;
Accession No.	Uniprot:P06748
Uniprot	P06748
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	33kDa
Storage	Store at 4°C in dark for 6 months

## Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

---

## Product Description

---

Involved in diverse cellular processes such as ribosome biogenesis, centrosome duplication, protein chaperoning, histone assembly, cell proliferation, and regulation of tumor suppressors p53/TP53 and ARF.

---

Note: This product is for in vitro research use only