

alpha + beta Synuclein Conjugated Antibody

Catalog No: #C56061



Package Size: #C56061-AF350 100ul #C56061-AF405 100ul #C56061-AF488 100ul
 #C56061-AF555 100ul #C56061-AF594 100ul #C56061-AF647 100ul
 #C56061-AF680 100ul #C56061-AF750 100ul #C56061-Biotin 100ul

Orders: order@signalwayantibody.com
 Support: tech@signalwayantibody.com

Description

Product Name	alpha + beta Synuclein Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Human Mouse Rat
Specificity	alpha + beta Synuclein Antibody detects endogenous levels of total alpha + beta Synuclein
Immunogen Description	A synthesized peptide derived from human alpha + beta Synuclein
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	alpha synuclein; beta synuclein ; NACP; PD1 ; SNCA; SNCB;
Accession No.	Uniprot:P37840/Q16143
Uniprot	P37840/Q16143
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	18kDa
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Product Description

Alpha-synuclein is a member of the synuclein family, which also includes beta- and gamma-synuclein. Synucleins are abundantly expressed in the brain and alpha- and beta-synuclein inhibit phospholipase D2 selectively.

Note: This product is for in vitro research use only