

HNF1 beta Conjugated Antibody

Catalog No: #C56108



Package Size: #C56108-AF350 100ul #C56108-AF405 100ul #C56108-AF488 100ul
 #C56108-AF555 100ul #C56108-AF594 100ul #C56108-AF647 100ul
 #C56108-AF680 100ul #C56108-AF750 100ul #C56108-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	HNF1 beta Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Human Mouse Rat
Specificity	HNF1 beta Antibody detects endogenous levels of HNF1 beta
Immunogen Description	A synthesized peptide derived from human HNF1 beta
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	Hepatocyte nuclear factor 1 beta; HNF 1B; HNF 2; HNF-1-beta; HNF-1B; HNF1 beta; HNF1 homeobox B; HNF1B; HNF2; HPC11; LF B3; MODY 5; TCF 2 protein; VHNF 1;
Accession No.	Uniprot:P35680
Uniprot	P35680
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	61kDa
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250
 AF405 conjugated: most applications: 1: 50 - 1: 250
 AF488 conjugated: most applications: 1: 50 - 1: 250
 AF555 conjugated: most applications: 1: 50 - 1: 250
 AF594 conjugated: most applications: 1: 50 - 1: 250
 AF647 conjugated: most applications: 1: 50 - 1: 250
 AF680 conjugated: most applications: 1: 50 - 1: 250
 AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Product Description

Belongs to the HNF1 homeobox family. Transcription factor, probably binds to the inverted palindrome 5'-GTTAATNATTAAC-3'.

Note: This product is for in vitro research use only