

Fibulin 5 Conjugated Antibody

Catalog No: #C56191



Package Size: #C56191-AF350 100ul #C56191-AF405 100ul #C56191-AF488 100ul
 #C56191-AF555 100ul #C56191-AF594 100ul #C56191-AF647 100ul
 #C56191-AF680 100ul #C56191-AF750 100ul #C56191-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Fibulin 5 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Human Mouse
Specificity	Fibulin 5 Antibody detects endogenous levels of total Fibulin 5
Immunogen Description	A synthesized peptide derived from human Fibulin 5
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	Fibulin-5; FIBL-5; Fibulin 5; FIBL 5; Fibulin5; FIBL5; Dance; Urine p50 protein; UP50; DANCE;
Accession No.	Uniprot:Q9UBX5
Uniprot	Q9UBX5
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	50kDa
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Product Description

Promotes adhesion of endothelial cells through interaction of integrins and the RGD motif. Could be a vascular ligand for integrin receptors and may play a role in vascular development and remodeling.

Note: This product is for in vitro research use only