

PDE4D Conjugated Antibody

Catalog No: #C56238



Package Size: #C56238-AF350 100ul #C56238-AF405 100ul #C56238-AF488 100ul
 #C56238-AF555 100ul #C56238-AF594 100ul #C56238-AF647 100ul
 #C56238-AF680 100ul #C56238-AF750 100ul #C56238-Biotin 100ul

Orders: order@signalwayantibody.com
 Support: tech@signalwayantibody.com

Description

Product Name	PDE4D Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Human
Specificity	PDE4D Antibody detects endogenous levels of total PDE4D
Immunogen Description	A synthesized peptide derived from human PDE4D
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	cAMP-specific 3',5'-cyclic phosphodiesterase 4D; DPDE3; PDE43; PDE4D; DPDE3;
Accession No.	Uniprot:Q08499
Uniprot	Q08499
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	91kDa
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Product Description

Hydrolyzes the second messenger cAMP, which is a key regulator of many important physiological processes.

Note: This product is for in vitro research use only