

# Trefoil Factor 3 Conjugated Antibody

Catalog No: #C56250



Package Size: #C56250-AF350 100ul #C56250-AF405 100ul #C56250-AF488 100ul  
 #C56250-AF555 100ul #C56250-AF594 100ul #C56250-AF647 100ul  
 #C56250-AF680 100ul #C56250-AF750 100ul #C56250-Biotin 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
 Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                       |  |
|-----------------------|--|
| Product Name          | Trefoil Factor 3 Conjugated Antibody   |
| Host Species          | Rabbit   |
| Clonality             | Monoclonal   |
| Isotype               | Rabbit IgG   |
| Purification          | Affinity-chromatography  |
| Species Reactivity    | Human  |
| Specificity           | Trefoil Factor 3 Antibody detects endogenous levels of total Trefoil Factor 3  |
| Immunogen Description | A synthesized peptide derived from human Trefoil Factor 3  |
| Conjugates            | Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750   |
| Other Names           | hITF; hP1.B; Intestinal trefoil factor; ITF; mITF; P1B; Polypeptide P1.B; TFF3; TFI; TREFOIL; Trefoil factor 3;  |
| Accession No.         | Uniprot:Q07654   |
| Uniprot               | Q07654   |
| Excitation Emission   | AF350: 346nm/442nm<br>AF405: 401nm/421nm<br>AF488: 493nm/519nm<br>AF555: 555nm/565nm<br>AF594: 591nm/614nm<br>AF647: 651nm/667nm<br>AF680: 679nm/702nm<br>AF750: 749nm/775nm |
| Calculated MW         | 9kDa   |
| Storage               | Store at 4°C in dark for 6 months  |

## Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

## Product Description

---

Involved in the maintenance and repair of the intestinal mucosa. Promotes the mobility of epithelial cells in healing processes (motogen).

---

Note: This product is for in vitro research use only