

## Thyroglobulin Conjugated Antibody

Catalog No: #C56331



Package Size: #C56331-AF350 100ul #C56331-AF405 100ul #C56331-AF488 100ul

#C56331-AF555 100ul #C56331-AF594 100ul #C56331-AF647 100ul

#C56331-AF680 100ul #C56331-AF750 100ul #C56331-Biotin 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Thyroglobulin Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Human Mouse Rat
Specificity	Thyroglobulin Antibody detects endogenous levels of total Thyroglobulin
Immunogen Description	A synthesized peptide derived from human Thyroglobulin
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	AITD3; hTG; TDH3; Tg; TGN; Thyroglobulin;
Accession No.	Uniprot:P01266
Uniprot	P01266
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	305kDa
Storage	Store at 4°C in dark for 6 months

## Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

## Product Description

---

Precursor of the iodinated thyroid hormones thyroxine (T4) and triiodothyronine (T3).

---

Note: This product is for in vitro research use only