

K Cadherin Conjugated Antibody

Catalog No: #C56356

Package Size: #C56356-AF350 100ul #C56356-AF405 100ul #C56356-AF488 100ul
 #C56356-AF555 100ul #C56356-AF594 100ul #C56356-AF647 100ul
 #C56356-AF680 100ul #C56356-AF750 100ul #C56356-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	K Cadherin Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Human Mouse Rat
Specificity	K Cadherin Antibody detects endogenous levels of total K Cadherin
Immunogen Description	A synthesized peptide derived from human K Cadherin
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	CAD6; Cdh6; K-cadherin; KCAD; Kidney cadherin;
Accession No.	Uniprot:P55285
Uniprot	P55285
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	130kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB:1:500~1:2000

IP:1:50

Product Description

Cadherins are calcium dependent cell adhesion proteins. They preferentially interact with themselves in a homophilic manner in connecting cells; cadherins may thus contribute to the sorting of heterogeneous cell types.

Note: This product is for in vitro research use only