

FSH beta $\alpha\Omega\frac{1}{2}\alpha\Omega\frac{1}{2}$ FSHBo $\Omega\frac{1}{2}\alpha\Omega\frac{1}{2}$ Conjugated Antibody

Catalog No: #C56398

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #C56398-AF350 100ul #C56398-AF405 100ul #C56398-AF488 100ul

#C56398-AF555 100ul #C56398-AF594 100ul #C56398-AF647 100ul

#C56398-AF680 100ul #C56398-AF750 100ul #C56398-Biotin 100ul

Description

Product Name	FSH beta $\alpha\Omega\frac{1}{2}\alpha\Omega\frac{1}{2}$ FSHBo $\Omega\frac{1}{2}\alpha\Omega\frac{1}{2}$ Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Human Rabbit
Specificity	FSH beta $\alpha\Omega\frac{1}{2}\alpha\Omega\frac{1}{2}$ FSHBo $\Omega\frac{1}{2}\alpha\Omega\frac{1}{2}$ Antibody detects endogenous levels of total FSH beta $\alpha\Omega\frac{1}{2}\alpha\Omega\frac{1}{2}$ FSHBo $\Omega\frac{1}{2}\alpha\Omega\frac{1}{2}$
Immunogen Description	A synthesized peptide derived from human FSH beta $\alpha\Omega\frac{1}{2}\alpha\Omega\frac{1}{2}$ FSHBo $\Omega\frac{1}{2}\alpha\Omega\frac{1}{2}$
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	Follicle Stimulating Hormon; Follitropin beta chain; Follitropin subunit beta; FSHB; FSHbeta;
Accession No.	Uniprot:P01225
Uniprot	P01225
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	22kDa
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Product Description

Follicle stimulating hormone enables ovarian folliculogenesis to the antral follicle stage and is essential for Sertoli cell proliferation and maintenance of sperm quality in the testis.

Note: This product is for in vitro research use only