

## 5 Lipoxygenase Conjugated Antibody

Catalog No: #C56414



Package Size: #C56414-AF350 100ul #C56414-AF405 100ul #C56414-AF488 100ul  
 #C56414-AF555 100ul #C56414-AF594 100ul #C56414-AF647 100ul  
 #C56414-AF680 100ul #C56414-AF750 100ul #C56414-Biotin 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
 Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	5 Lipoxygenase Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Human Mouse Rat
Specificity	5 Lipoxygenase Antibody detects endogenous levels of total 5 Lipoxygenase
Immunogen Description	A synthesized peptide derived from human 5 Lipoxygenase
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	ALOX5; 5-LO; 5-LOX; 5LPG; LOG5; MGC163204; LOX5; 5-lipoxygenase ; 5 Lipoxygenase; 5 LOX; ALOX 5;
Accession No.	Uniprot:P09917
Uniprot	P09917
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	78kDa
Storage	Store at 4°C in dark for 6 months

## Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

## Product Description

---

Catalyzes the first step in leukotriene biosynthesis, and thereby plays a role in inflammatory processes.

---

Note: This product is for in vitro research use only