

Human Kappa Chain Conjugated Antibody

Catalog No: #C56434



Package Size: #C56434-AF350 100ul #C56434-AF405 100ul #C56434-AF488 100ul
 #C56434-AF555 100ul #C56434-AF594 100ul #C56434-AF647 100ul
 #C56434-AF680 100ul #C56434-AF750 100ul #C56434-Biotin 100ul

Orders: order@signalwayantibody.com
 Support: tech@signalwayantibody.com

Description

| | |
|-----------------------|--|
| Product Name | Human Kappa Chain Conjugated Antibody |
| Host Species | Rabbit |
| Clonality | Monoclonal |
| Isotype | Rabbit IgG |
| Purification | Affinity-chromatography |
| Species Reactivity | Human |
| Specificity | Human Kappa Chain Antibody detects endogenous levels of total Human Kappa Chain |
| Immunogen Description | A synthesized peptide derived from human Human Kappa Chain |
| Conjugates | Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750 |
| Other Names | HCAK1; Ig kappa chain C region; IGKC; IMMUNOGLOBULIN InV; Km; |
| Accession No. | Uniprot:P01834 |
| Uniprot | P01834 |
| Excitation Emission | AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm |
| Calculated MW | 26kDa |
| Storage | Store at 4°C in dark for 6 months |

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Product Description

Immunoglobulins belong to a group of related glyco proteins which make up 20% of serum proteins. Antigens and immunoglobulins react to confer immunity to individuals. Immunoglobulins have similar structures of two identical heavy chains and two identical light chains.

Note: This product is for in vitro research use only