

Glutathione Peroxidase 2 Conjugated Antibody

Catalog No: #C56437



Package Size: #C56437-AF350 100ul #C56437-AF405 100ul #C56437-AF488 100ul
 #C56437-AF555 100ul #C56437-AF594 100ul #C56437-AF647 100ul
 #C56437-AF680 100ul #C56437-AF750 100ul #C56437-Biotin 100ul

Orders: order@signalwayantibody.com
 Support: tech@signalwayantibody.com

Description

Product Name	Glutathione Peroxidase 2 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Human
Specificity	Glutathione Peroxidase 2 Antibody detects endogenous levels of total Glutathione Peroxidase 2
Immunogen Description	A synthesized peptide derived from human Glutathione Peroxidase 2
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	GI GPx; GPRP; GPRP-2; GPx-GI; GPX2; GSHPx 2; GSHPx GI; GSHPx-2; GSHPx-GI;
Accession No.	Uniprot:P18283
Uniprot	P18283
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	22kDa
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Product Description

Could play a major role in protecting mammals from the toxicity of ingested organic hydroperoxides. Tert-butyl hydroperoxide, cumene hydroperoxide and linoleic acid hydroperoxide but not phosphatidylcholine hydroperoxide, can act as acceptors.

Note: This product is for in vitro research use only