

## TNFSF9 Conjugated Antibody

Catalog No: #C56562



Package Size: #C56562-AF350 100ul #C56562-AF405 100ul #C56562-AF488 100ul  
 #C56562-AF555 100ul #C56562-AF594 100ul #C56562-AF647 100ul  
 #C56562-AF680 100ul #C56562-AF750 100ul #C56562-Biotin 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
 Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	TNFSF9 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Human
Specificity	TNFSF9 Antibody detects endogenous levels of total TNFSF9
Immunogen Description	A synthesized peptide derived from human TNFSF9
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	4 1BB L; 4 1BBL; 4-1BB ligand; 4-1BBL; Cd137l; Cd157l; Ly63l; Tnfsf9; TNLG5A;
Accession No.	Uniprot:P41273
Uniprot	P41273
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	25-32kDa
Storage	Store at 4°C in dark for 6 months

## Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

## Product Description

---

Cytokine that binds to TNFRSF9. Induces the proliferation of activated peripheral blood T-cells. May have a role in activation-induced cell death (AICD). May play a role in cognate interactions between T-cells and B-cells/macrophages.

---

Note: This product is for in vitro research use only