

EBP1 / PA2G4 Conjugated Antibody

Catalog No: #C56650



Package Size: #C56650-AF350 100ul #C56650-AF405 100ul #C56650-AF488 100ul
 #C56650-AF555 100ul #C56650-AF594 100ul #C56650-AF647 100ul
 #C56650-AF680 100ul #C56650-AF750 100ul #C56650-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	EBP1 / PA2G4 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Human Mouse Rat
Specificity	EBP1 / PA2G4 Antibody detects endogenous levels of total EBP1 / PA2G4
Immunogen Description	A synthesized peptide derived from human EBP1 / PA2G4
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	hG4 1; Mpp1; p38 2G4; Pa2g4; Plfap;
Accession No.	Uniprot:Q9UQ80
Uniprot	Q9UQ80
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	48kDa
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Product Description

May play a role in a ERBB3-regulated signal transduction pathway. Seems be involved in growth regulation. Acts a corepressor of the androgen receptor (AR) and is regulated by the ERBB3 ligand neuregulin-1/herregulin (HRG).

Note: This product is for in vitro research use only