

GLUD1 Conjugated Antibody

Catalog No: #C56733



Package Size: #C56733-AF350 100ul #C56733-AF405 100ul #C56733-AF488 100ul

#C56733-AF555 100ul #C56733-AF594 100ul #C56733-AF647 100ul

#C56733-AF680 100ul #C56733-AF750 100ul #C56733-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	GLUD1 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Human Mouse Rat
Specificity	GLUD1 Antibody detects endogenous levels of total GLUD1
Immunogen Description	A synthesized peptide derived from human GLUD1
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	GDH; Gdh-X; GDH1; GLUD; Glud1; Glud1a; Glud1b; Gludl; MRG-2;
Accession No.	Uniprot:P00367
Uniprot	P00367
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	50kDa
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Product Description

May be involved in learning and memory reactions by increasing the turnover of the excitatory neurotransmitter glutamate.

Note: This product is for in vitro research use only