

GCYA3 Conjugated Antibody

Catalog No: #C56889



Package Size: #C56889-AF350 100ul #C56889-AF405 100ul #C56889-AF488 100ul
 #C56889-AF555 100ul #C56889-AF594 100ul #C56889-AF647 100ul
 #C56889-AF680 100ul #C56889-AF750 100ul #C56889-Biotin 100ul

Orders: order@signalwayantibody.com
 Support: tech@signalwayantibody.com

Description

Product Name	GCYA3 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Human
Specificity	GCYA3 Antibody detects endogenous levels of total GCYA3
Immunogen Description	A synthesized peptide derived from human GCYA3
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	GUCY1A3; GC-S-alpha-1; GUC1A3; GUCA3; GCS-alpha-3; GUCY1A1; GC-SA3; GCS-alpha-1; GUCSA3;
Accession No.	Uniprot:Q02108
Uniprot	Q02108
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	73kDa
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Product Description

Belongs to the adenylyl cyclase class-4/guanylyl cyclase family. Contains 1 guanylate cyclase domain.

Note: This product is for in vitro research use only