

CC2D1A Conjugated Antibody

Catalog No: #C56901



Package Size: #C56901-AF350 100ul #C56901-AF405 100ul #C56901-AF488 100ul
 #C56901-AF555 100ul #C56901-AF594 100ul #C56901-AF647 100ul
 #C56901-AF680 100ul #C56901-AF750 100ul #C56901-Biotin 100ul

Orders: order@signalwayantibody.com
 Support: tech@signalwayantibody.com

Description

Product Name	CC2D1A Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Human Mouse Rat
Specificity	CC2D1A Antibody detects endogenous levels of total CC2D1A
Immunogen Description	A synthesized peptide derived from human CC2D1A
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	Coiled coil and C2 domain containing 1A; FREUD 1; Freud 1/Aki1; MRT3;
Accession No.	Uniprot:Q6P1N0
Uniprot	Q6P1N0
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	130kDa
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Product Description

May play a role in the altered regulation of HTR1A associated with anxiety and major depression. It also mediates HDAC-independent repression of HTR1A promoter in neuronal cell.

Note: This product is for in vitro research use only