

HTF9C Conjugated Antibody

Catalog No: #C56915



Package Size: #C56915-AF350 100ul #C56915-AF405 100ul #C56915-AF488 100ul
 #C56915-AF555 100ul #C56915-AF594 100ul #C56915-AF647 100ul
 #C56915-AF680 100ul #C56915-AF750 100ul #C56915-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	HTF9C Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Human
Specificity	HTF9C Antibody detects endogenous levels of total HTF9C
Immunogen Description	A synthesized peptide derived from HTF9C
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	HTF9C; TRMT2A;
Accession No.	Uniprot:Q8IZ69
Uniprot	Q8IZ69
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	75kDa
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Product Description

May be involved in nucleic acid metabolism and/or modifications.

Note: This product is for in vitro research use only