

Kallikrein 8 Conjugated Antibody

Catalog No: #C57027



Package Size: #C57027-AF350 100ul #C57027-AF405 100ul #C57027-AF488 100ul
 #C57027-AF555 100ul #C57027-AF594 100ul #C57027-AF647 100ul
 #C57027-AF680 100ul #C57027-AF750 100ul #C57027-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Kallikrein 8 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Human
Specificity	Kallikrein 8 Antibody detects endogenous levels of total Kallikrein 8
Immunogen Description	A synthesized peptide derived from human Kallikrein 8
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	BSP1; hK8; HNP; HPN; Kallikrein8; Kik8; NRPN; PRSS19; TADG14; UNQ283;
Accession No.	Uniprot:O60259
Uniprot	O60259
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	35kDa
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Product Description

Serine protease which is capable of degrading a number of proteins such as casein, fibrinogen, kininogen, fibronectin and collagen type IV. Also cleaves L1CAM in response to increased neural activity.

Note: This product is for in vitro research use only