

Calcipressin 1 Conjugated Antibody

Catalog No: #C57033



Package Size: #C57033-AF350 100ul #C57033-AF405 100ul #C57033-AF488 100ul
 #C57033-AF555 100ul #C57033-AF594 100ul #C57033-AF647 100ul
 #C57033-AF680 100ul #C57033-AF750 100ul #C57033-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Calcipressin 1 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Human Mouse Rat
Specificity	Calcipressin 1 Antibody detects endogenous levels of total Calcipressin 1
Immunogen Description	A synthesized peptide derived from human Calcipressin 1
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	Adapt78; Calcipressin 1; CSP1; DSC1; DSCR1; MCIP1; RCAN1; RCN1;
Accession No.	Uniprot:P53805
Uniprot	P53805
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	28kDa
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Product Description

Inhibits calcineurin-dependent transcriptional responses by binding to the catalytic domain of calcineurin A. Could play a role during central nervous system development.

Note: This product is for in vitro research use only