

ICAM3 / CD50 Conjugated Antibody

Catalog No: #C57053



Package Size: #C57053-AF350 100ul #C57053-AF405 100ul #C57053-AF488 100ul
 #C57053-AF555 100ul #C57053-AF594 100ul #C57053-AF647 100ul
 #C57053-AF680 100ul #C57053-AF750 100ul #C57053-Biotin 100ul

Orders: order@signalwayantibody.com
 Support: tech@signalwayantibody.com

Description

Product Name	ICAM3 / CD50 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Human
Specificity	ICAM3 / CD50 Antibody detects endogenous levels of total ICAM3 / CD50
Immunogen Description	A synthesized peptide derived from human ICAM3 / CD50
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	CD50; CDw50; ICAM 3; ICAM3; ICAMR; Intercellular adhesion molecule 3;
Accession No.	Uniprot:P32942
Uniprot	P32942
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	100-160kDa
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Product Description

ICAM proteins are ligands for the leukocyte adhesion protein LFA-1 (integrin alpha-L/beta-2). ICAM3 is also a ligand for integrin alpha-D/beta-2.

Note: This product is for in vitro research use only