

CD32A + CD32B + CD32C Conjugated Antibody

Catalog No: #C57085



Package Size: #C57085-AF350 100ul #C57085-AF405 100ul #C57085-AF488 100ul
 #C57085-AF555 100ul #C57085-AF594 100ul #C57085-AF647 100ul
 #C57085-AF680 100ul #C57085-AF750 100ul #C57085-Biotin 100ul

Orders: order@signalwayantibody.com
 Support: tech@signalwayantibody.com

Description

Product Name	CD32A + CD32B + CD32C Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Human
Specificity	CD32A + CD32B + CD32C Antibody detects endogenous levels of total CD32A + CD32B + CD32C
Immunogen Description	A synthesized peptide derived from human CD32A + CD32B + CD32C
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	CD32; CD32A; CD32B; CD32C; CDw32; FCG2; FcGR; FCGR2; FCGR2A; FCGR2A1; FCGR2B; FCGR2C; Fcgr3; Fcgr3a; FcgRII; FcRII a; FcRII b; FcRII c; FcRIIC; IGFR2;
Accession No.	Uniprot:P12318/P31994/P31995
Uniprot	P12318/P31994/P31995
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	37-43kDa
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250
 AF405 conjugated: most applications: 1: 50 - 1: 250
 AF488 conjugated: most applications: 1: 50 - 1: 250
 AF555 conjugated: most applications: 1: 50 - 1: 250
 AF594 conjugated: most applications: 1: 50 - 1: 250
 AF647 conjugated: most applications: 1: 50 - 1: 250
 AF680 conjugated: most applications: 1: 50 - 1: 250
 AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Product Description

Binds to the Fc region of immunoglobulins gamma. Low affinity receptor. By binding to IgG it initiates cellular responses against pathogens and soluble antigens. Promotes phagocytosis of opsonized antigens.

Note: This product is for in vitro research use only