

Anti-CD8a Antibody

Catalog No: #33199



Package Size: #33199-1 50ul #33199-2 50ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Anti-CD8a Antibody
Host Species	Rat
Clonality	Monoclonal
Purification	The antibody was purified by affinity chromatography.
Applications	IF, FC
Species Reactivity	Mouse
Immunogen Description	Mouse thymus or spleen
Other Names	MAL; T-cell surface glycoprotein CD8 alpha chain; T-lymphocyte differentiation antigen T8/Leu-2; CD8a
Accession No.	Swiss Prot:P01731GeneID:12525
Uniprot	P01731
GeneID	12525
Formulation	Rat IgG2a kappa. Liquid in PBS, pH 7.3, and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

IF:1:50-1:200

Note: This product is for in vitro research use only